



Euro Disposables
C/ Jose Ortega y Gasset 13, 7^a, 28006 Madrid, SPAIN
Tel: +34 609 741518, e-mail: info@eurodisposables.com

Disposable Auto-disable Syringe

General features

Size: 0.5ml, 1ml, 2.5ml, 3ml, 5ml, 10ml

Sterile, non-toxic, non-pyrogen, single use only

Clear scale mark, quality stainless steel needle

Prevent the re-use of syringes and thus the transmission of diseases

Reduces waste of medication to a minimum

Shelf life: 5years

Quality certificates: CE/ISO9001/ISO13485

We also supply OEM services

Characteristics of the components

Syringe structure: Barrel, Plunger, Piston, Needle, Nozzle cap, Needle file, Stop block and The upper connecting base

Material:

Barrel: made of high transparent medical grade PP

Plunger: made of high transparent medical grade PP

Piston: made of natural rubber with two retaining rings or synthetic non-cytotoxic rubber

Needle: high quality stainless steel SUS304

Nozzle cap: made of high transparent medical grade PP

Needle file: made of high transparent medical grade PP

Stop block: made of high transparent medical grade PP

The upper connecting base: made of high transparent medical grade PP

Graduation : indelible ink

Blister or plastic package: medical grade paper and plastic film

Individually packing: PE bag or blister packaging

Inner packing: box/polybag

Outer packing: carton

Sterilization method: EO gas

