



Euro Disposables- Import/Export International Business & Market Development e.U.

EU Identification Number: ATU68505129

Dk.-Theodor- Gugl- Weg 16, 5081 Anif AUSTRIA

ONE USE Medical Products, CE Labelled

1.1 Disposables Infusion Sets

1.2 Urine Bags

1.3 Disposables Feeding Tube

1.4 Disposables Stomach

1.5 Rectal Tube Nelaton catheter

1.6 Disposable Auto-disable Syringe

1.7 Disposables Sterile Syringes- Luer Slip

1.8 Disposables Syringes 2-Parts

1.9 Disposables Syringes 3-Parts

Other products please request specifications.



Disposable Infusion Set

➤ **General features**

Sterile, Non-Pyrogenic, Non-toxic

Used for body vein infusion

Shelf life: 3years

Quality certificates: CE/ISO9001/ISO13485

We also supply OEM services

➤ **Characteristics of the components**

Spike protector: made of PE, air-vent or Non-vent

Spike : made of ABS, Sharp and easy

Drip chamber: made of soft PVC, Blowing mold

Flow regulator: made of HDPE

Tube: made of PVC

Male conical fitting: made of HDPE

“Y” injection port: medical latex or latex free,

With or without injection needle


Blister or plastic package: medical grade paper and plastic film

Individually packing: PE bag or blister packaging

Inner packing: box/polybag

Outer packing: carton

➤ **Sterilization method: EO gas**





Urine bag

➤ General features

Size: 2000ml, 1500ml, 1000ml etc (we also can produce the urine bag according to your request)

Non-toxic, for single use, mainly use for liquid-leading and urine collection after operation

Quality certificates: CE/ISO9001/ISO13485

We also supply OEM services

➤ Characteristics of the components

Made of medical grade PVC material

90cm long super smooth kink resistant tubing is easy to milk and ensures unrestricted flow

Non return valve prevents back slow of bacteria migration

Conic trapped connector


With or without pull-push valve, T valve, screw valve

Plastic package: medical plastic film

Individually packing: PE bag

Outer packing: carton

➤ Sterilization method: EO gas





Euro Disposables- Import/Export International Business & Market Development e.U.

EU Identification Number: ATU68505129

Dk.-Theodor- Gugl- Weg 16, 5081 Anif AUSTRIA

Disposable feeding tube

General features

Size: F4, F8, F10, F12, F14, F16, F18, F20, F22

Transparent, disposable, non-toxic, single package

Through the patient's nostril to irrigate nutrition, without cross contamination

Length marking facilitate to insertion

Quality certificates: CE/ISO9001/ISO13485

We also supply OEM services

Characteristics of the components

Made of non-toxic PVC, medical grade

Two lateral eyes & open tip

Child: 40cm adult: 120cm

Funnel connector with cap

Contain DEHP

Blister or plastic package: medical grade paper and plastic film

Individually packing: PE bag or blister packaging

Inner packing: box/polybag

Outer packing: carton

Sterilization method: EO gas





Euro Disposables- Import/Export International Business & Market Development e.U.

EU Identification Number: ATU68505129

Dk.-Theodor- Gugl- Weg 16, 5081 Anif AUSTRIA

Disposables Tube for the Stomach

General features

Size: F6, F8, F10, F12, F14, F16, F18, F20, F22

To help with nutriment injecting, phlegm and residual conducting to patient who cannot take food

Temporarily after surgical operation

Normative color coded connector, easy to recognize

Quality certificates: CE/ISO9001/ISO13485

We also supply OEM services

Characteristics of the components

Made of non-toxic PVC, medical grade

Four lateral eyes & open tip

Length : 120cm

Funnel connector with X-radiopaque graduation

Contain DEHP

Blister or plastic package: medical grade paper and plastic film

Individually packing: PE bag or blister packaging

Inner packing: box/polybag

Outer packing: carton

Sterilization method: EO gas





Euro Disposables- Import/Export International Business & Market Development e.U.

EU Identification Number: ATU68505129

Dk.-Theodor- Gugl- Weg 16, 5081 Anif AUSTRIA

Rectal Tube

General features

Size: F24, F26, F28, F30, F32, F34, F36

For introduction of enem a solution into rectum to release

Soft, frosted and kink resistant tubing

Proximal end is fitted with Universal Funnel Shaped connector for sale connection to Urine bag

Colour coded plain connector for easy identification of size

Quality certificates: CE/ISO9001/ISO13485

We also supply OEM services

Characteristics of the components

Made of non-toxic PVC, medical grade

Two lateral eyes & rounded closed tip

Transparent & smooth

Length : 40cm

Funnel connector

Contain DEHP

Blister or plastic package: medical grade paper and plastic film

Individually packing: PE bag or blister packaging

Inner packing: box/polybag

Outer packing: carton

Sterilization method: EO gas





Euro Disposables- Import/Export International Business & Market Development e.U.

EU Identification Number: ATU68505129

Dk.-Theodor- Gugl- Weg 16, 5081 Anif AUSTRIA

Nelaton catheter

General features

Size: F6, F8, F10, F12, F14, F16, F18, F20, F22

Soft, frosted and kink resistant tubing

Proximal end is fitted with Universal Funnel Shaped connector for safe connection to Urine bag

For short term baladder catheterisation through urethra in males

Colour coded plain connector for easy identification of size

Quality certificates: CE/ISO9001/ISO13485

We also supply OEM services

Characteristics of the components

Made of non-toxic PVC, medical grade

Two lateral eyes & rounded closed tip

Transparent & smooth

Length : 40cm

Funnel connector

Contain DEHP


Blister or plastic package: medical grade paper and plastic film

Individually packing: PE bag or blister packaging

Inner packing: box/polybag

Outer packing: carton

Sterilization method: EO gas





Disposable Auto-disable Syringe

➤ **General features**

Size: 0.5ml, 1ml, 2.5ml, 3ml, 5ml, 10ml

Sterile, non-toxic, non-pyrogen, single use only

Clear scale mark, quality stainless steel needle

Prevent the re-use of syringes and thus the transmission of diseases

Reduces waste of medication to a minimum

Shelf life: 5years

Quality certificates: CE/ISO9001/ISO13485

We also supply OEM services

➤ **Characteristics of the components**

Syringe structure: Barrel, Plunger, Piston, Needle, Nozzle cap, Needle file, Stop block and The upper connecting base

➤ **Material:**

Barrel: made of high transparent medical grade PP

Plunger: made of high transparent medical grade PP

Piston: made of natural rubber with two retaining rings or synthetic non-cytotoxic rubber

Needle: high quality stainless steel SUS304

Nozzle cap: made of high transparent medical grade PP

Needle file: made of high transparent medical grade PP

Stop block: made of high transparent medical grade PP

The upper connecting base: made of high transparent medical grade PP

Graduation : indelible ink

Blister or plastic package: medical grade paper and plastic film

Individually packing: PE bag or blister packaging

Inner packing: box/polybag

Outer packing: carton

➤ **Sterilization method:** EO gas





Disposables Sterile Syringes- Luer Slip

- General features

Injection-tube size: 1ml,2ml/2.5ml,3ml,5ml,10ml,20ml,30ml,50/60ml

Sterile, non-toxic, non-pyrogenic, single use.

3 parts, Luer lock

Material for barrel and plunger: medical and high transparent PP.

With or without needle.

Hypodermic needle:29G,28G,27G,26G,25G,24G,23G,22G,21G,20G,18G,16G

Material for gasket: Natural latex / Synthetic rubber(latex free).

Transparent casing for easy observance; strongly adhesive printing ink of the scale that will not fall off.

Certificates: ISO&CE.

We also supply OEM services

- Characteristics of the components

Barrel made of high transparent medical grade PP

Plunger made of high transparent medical grade PP

Gasket made of nature rubber latex or synthetic rubber

Graduation: indelible ink

Blister or plastic package: medical grade paper and plastic film

Individually packing: PE bag or blister packaging

Inner packing: box/polybag

Outer packing: carton

- Sterilization method: EO gas





Euro Disposables- Import/Export International Business & Market Development e.U.

EU Identification Number: ATU68505129

Dk.-Theodor- Gugl- Weg 16, 5081 Anif AUSTRIA

Disposables Syringes 2-Parts

General features

Injection-tube size: 2ml, 5ml, 10ml, 20ml

Sterile, non-toxic, non-pyrogenic, for single use.

2-parts luer slip

Material for barrel and plunger: medical and high transparent PP.

With or without needle.

Hypodermic needle : 16G, 18G, 20G, 21G, 22G, 23G, 24G, 25G, 26G, 27G, 28G, 29G

Transparent casing for easy observance; strongly adhesive printing ink of the scale that will not fall off.

Certificates: ISO&CE

We also supply OEM services

Characteristics of the components

Barrel made of high transparent medical grade PP

Plunger made of high transparent medical grade PP

Graduation: indelible ink

Blister or plastic package: medical grade paper and plastic film

Individually packing: PE bag or blister packaging

Inner packing: box/polybag

Outer packing: carton

Sterilization method: EO gas





Disposables Syringes 3-Parts

➤ **General features**

Injection-tube size:1ml 2ml 3ml 5ml 10ml 20ml 30ml 50ml 60ml

Sterile, non-toxic, non-pyrogenic, for single use.

3 parts,Luer slip or Luer lock .

Material for barrel and plunger: medical and high transparent PP.

With or without needle.

Hypodermic needle:29G,28G,27G,26G,25G,24G,23G,22G,21G,20G,18G,16G

Material for gasket: Natural latex / Synthetic rubber(latex free).

Transparent casing for easy observance; strongly adhesive printing ink of the scale that will not fall off.

Certificates: ISO&CE.

We also supply OEM services

➤ **Characteristics of the components**

Barrel made of high transparent medical grade PP

Plunger made of high transparent medical grade PP

Gasket made of nature rubber latex or synthetic rubber

Graduation: indelible ink

Blister or plastic package: medical grade paper and plastic film

Individually packing: PE bag or blister packaging

Inner packing: box/polybag

Outer packing: carton

➤ **Sterilization method: EO gas**

