


Infusion set REF:JSM1004

- ABS spike without air-vent,
- Big drip chamber 62mm
- Filter (15 μ)
- 20 Drop/ml
- Safety roll clamp
- Latex free bulb
- Luer lock connector with cap
- Tube 150cm
- Paper peel package
- 25pcs/Bag,100pcs/carton,
- Carton size:41x23.5x29.5cm

- Pyrogen free
- Latex free
- For single use
- ETO sterilized



| N° | Name | Material |
|----|-------------------|-----------------|
| 1 | Protective cap | PE |
| 2 | Vented spike | ABS |
| 3 | Drip chamber | PVC |
| 4 | Fluid filter | ABS & Nylon |
| 5 | Roller clamp | ABS |
| 6 | Latex free | Latex free |
| 7 | Male LL with cap | Transparent ABS |
| 8 | Tube 3x4.1x1500mm | PVC |