



# ISOLATION EYES COVER





**[Model / Specification] GLYZ1-1**

---

## PRODUCT DESCRIPTION

---

It consists of a protective cover,  
Nonsterile,  
single-use.



PRODUCT DISPLAY



PRODUCT APPLICATION

For personal face protection.



PRODUCT DISPLAY





## SPECIFICATIONS



Original: China


Terms of trade: TT

1pcs product size	203.5*57.5*86mm
1pcs product weight	51.5g
package	The Chinese and English instructions are printed on the package
1 boxes quantity	20pcs
carton quantity	200pcs
N.W	10.3kg
G.W	16.52kg
Meast	950*480*375mm
Certificate	CE、CFDA



### INSTRUCTIONS

1. Unpack the inner packaging, remove the isolation mask product, and remove the protective film on both sides.
2. After the two ends of the foam strip (commonly known as elastic band) are fastened, wear it on the forehead to adjust it comfortably.



### ATTENTION

1. Do not expose the isolation eyes cover to strong sunlight.
2. Handle carefully. The part of the isolation eyes cover must not touch sharp objects.