

**WARNING**

Do not use if the packaging damaged or time expires.

**STORAGE CONDITIONS**

It should be stored in a room with relative humidity not exceeding 80%, no corrosiveas, and well-ventilated.

# ISOLATION FACE MASK

**[Model ]GLMZ1-1**



## PRODUCT DESCRIPTION

It consists of a protective cover,  
foam strip and fixing  
device made of  
polymer materials.  
Nonsterile,  
single-use.



PRODUCT DISPLAY



PRODUCT DISPLAY



## SPECIFICATIONS



Original: China


Terms of trade: TT

1pcs product size	350*240*30mm
1pcs product weight	60.7g
package	The Chinese and English instructions are printed on the package
1 boxes quantity	No inner box
carton quantity	200pcs
N.W	12.14kg
G.W	13.87kg
Meast	680*383*415mm
Certificate	CE、CFDA



## INSTRUCTIONS

- 1、Unpack the inner packaging, remove the isolation mask product, and remove the protective film on both sides.
- 2、After the two ends of the foam strip (commonly known as elastic band) are fastened, wear it on the forehead to adjust it comfortably.



## ATTENTION

- 1、Do not expose the isolation face mask to strong sunlight.
- 2、Handle with care. The part of the isolation face mask must not touch sharp objects.